



PROBE N°19

Ref. IDSOND19



Exemples de SONDES RHEIN – Holtex
Réf. IDSOND01, IDSOND19, IDSON23
Existe également en n°2, 6, 8, 9, 16 et 17



19/06/2026 09:04:52 – For information purposes only, can be subject to change

PROBE N°19

Ref. IDSOND19



Description

Rhein probes

9 different models are available, contact us

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104738157	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	0.9
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.9
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16.2
Packed unit / Weight	0.016	Bare Unit / Weight (kg)	0.016

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---