

PROBE N°23

Ref. IDSOND23



CE



Description

Rhein probes

9 different models are available, contact us

19/06/2026 08:30:02 – For information purposes only, can be subject to change

PROBE N°23

Ref. IDSOND23



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104738164	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	14.7
Packed unit / Weight	0.023	Bare Unit / Weight (kg)	0.023
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 19/06/2026 08:30:02 – For information purposes only, can be subject to change