

# HARTMANN RETRACTOR STAINLESS STEEL WIRE 16 CM (PAIR)



Ref. IE03616

## Features

CE Marking	CE	Medical device	Class I
Bar code	08961104721715	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	5.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.1	Bare Unit / Weight (kg)	0.085
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 18/06/2026 22:01:59 – For information purposes only, can be subject to change