

LOVE RETRACTOR 45° ANGULATED, 19CM

Ref. IE05401



Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104721883	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	4.8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	18.2
Packed unit / Weight	0.025	Bare Unit / Weight (kg)	0.025
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 19/06/2026 01:22:50 – For information purposes only, can be subject to change