

# ALM RETRACTOR 10 CM

Ref. IE11110



## Features

|                                    |                |                         |          |
|------------------------------------|----------------|-------------------------|----------|
| CE Marking                         | CE             | Medical device          | Class I  |
| Bar code                           | 08961104722835 | HS Code                 | 90189084 |
| Warranty (months)                  | 24             | Country                 | Pakistan |
| <b>DIMENSIONS AND WEIGHT</b>       |                |                         |          |
| Packed unit / Width (cm)           | nc             | Bare Unit / Width (cm)  | 3.1      |
| Packed unit / Height (cm)          | nc             | Bare Unit / Height (cm) | 0.7      |
| Packed unit / Length (cm)          | nc             | Bare Unit / Length (cm) | 10       |
| Packed unit / Weight               | 0.019          | Bare Unit / Weight (kg) | 0.019    |
| <b>SPECIFIC STORAGE CONDITIONS</b> |                |                         |          |
| Minimum temperature                | NA             | Max temperature         | NA       |
| Minimum humidity                   | NA             | Maximum humidity        | NA       |
| <b>OTHER FEATURES</b>              |                |                         |          |
| Packed unit / Quantity             | 1              |                         |          |

📅 18/06/2026 22:19:21 – For information purposes only, can be subject to change