

BECKMANN RETRACTOR (4X4 TEETH) 30 CM



Ref. IE13430

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104722965	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	12.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	30.5
Packed unit / Weight	0.22	Bare Unit / Weight (kg)	0.22
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 19/06/2026 01:05:20 – For information purposes only, can be subject to change