

TUFFIER RETRACTOR 50 X 40 MM**Ref. IE29540****Features**

CE Marking	CE	Medical device	Class I
Bar code	08961104731578	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	12.6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	23.8
Packed unit / Weight	0.465	Bare Unit / Weight (kg)	0.462

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---