

YASARGIL/LEYLA SELF-RETAINING RETRACTOR FOR SKULL MOUNTING

Ref. IE32200



19/06/2026 00:19:12 – For information purposes only, can be subject to change

YASARGIL/LEYLA SELF-RETAINING RETRACTOR FOR SKULL MOUNTING

Ref. IE32200



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104726222	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.45
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		