

ROCHESTER-PEAN FORCEPS STRAIGHT NO TEETH 15 CM

Ref. IP02215



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL

19/06/2026 01:19:59 – For information purposes only, can be subject to change

ROCHESTER-PEAN FORCEPS STRAIGHT NO TEETH 15 CM



Ref. IP02215

Convenient to stop bleeding without damaging fragile tissues during an operation. Also used to perform care (dressings, etc.).

Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations.

Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016).

- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" leaflet

Features

| | | | |
|------------------------------------|----------------|-------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 08961104749641 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Bare Unit / Width (cm) | 7.5 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 0.7 |
| Packed unit / Length (cm) | nc | Bare Unit / Length (cm) | 15.6 |
| Packed unit / Weight | 0.044 | Bare Unit / Weight (kg) | 0.043 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

📅 19/06/2026 01:19:59 – For information purposes only, can be subject to change