

LONDON COLLEGE TWEEZERS 15 CM

Ref. IP04715



CE



London College Tweezers 15 cm

Description

London College Tweezers 15 cm

19/06/2026 02:15:52 – For information purposes only, can be subject to change

LONDON COLLEGE TWEEZERS 15 CM

Ref. IP04715



Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104714373	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.6
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15
Packed unit / Weight	0.023	Bare Unit / Weight (kg)	0.023
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 19/06/2026 02:15:52 – For information purposes only, can be subject to change