


BACKHAUS TOWEL FORCEPS 14 CM

Ref. IP10614



19/06/2026 01:21:24 – For information purposes only, can be subject to change

BACKHAUS TOWEL FORCEPS 14 CM

Ref. IP10614



To attach/secure surgical towels or other devices around the surgical incision.

Description

Holtex is specialist in stainless steel surgical instruments since 1980. Holtex instruments are all made of top quality stainless steel in order to guarantee good resistance to corrosion and durability.

> This reference is made with Stainless Steel AISI 410.

The Holtex surgical instruments are classified as medical devices and comply ISO 13485 standards and the 93/42/EC directive.

EN ISO 13485: 2016, ISO 14971: 2019, EN 1041: 2008, and EN ISO 15223-1

Holtex Instruments are guaranteed for 2 years under normal conditions of use. The warranty does not cover wear parts, shocks and maintenance defects made by the user

Packaging and traceability: Holtex instruments are supplied in individual bags, sealed and labelled with : Logo, and address of the manufacturer, the product reference, the detailed designation in French and English, the lot number ("LOT"), the CE marking, and the barcode.

Cleaning: Please read carefully instructions of cleaning and sterilisation of stainless-steel equipment inserted in each individual bag. Instructions can be also downloaded from the product page on our website.

BACKHAUS TOWEL FORCEPS 14 CM

Ref. IP10614



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104720367	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	14.2
Packed unit / Weight	0.03	Bare Unit / Weight (kg)	0.03
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 19/06/2026 01:21:24 – For information purposes only, can be subject to change