

## MURPHY-PEAN FORCEPS STRAIGHT 13 CM

Ref. IP21413



CE

DISPOSITIF MEDICAL  
CLASSE  
**I**  
MEDICAL DEVICE

GARANTIE  
2 ANS  
**2**  
2 YEARS  
GUARANTEE

STAINLESS  
**INOX**  
STEEL

18/06/2026 23:52:41 – For information purposes only, can be subject to change

# MURPHY-PEAN FORCEPS STRAIGHT 13 CM

Ref. [IP21413](#)



**Used during surgical procedures (stop bleeding, hold compress, etc.)**

## Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations.

- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" leaflet

## Features

|                   |                |                |          |
|-------------------|----------------|----------------|----------|
| CE Marking        | CE             | Medical device | Class I  |
| Bar code          | 08961104721050 | HS Code        | 90189084 |
| Warranty (months) | 24             | Country        | Pakistan |

### DIMENSIONS AND WEIGHT

|                           |      |                           |       |
|---------------------------|------|---------------------------|-------|
| Packed unit / Width (cm)  | nc   | Packed unit / Height (cm) | nc    |
| Packed unit / Length (cm) | nc   | Bare Unit / Length (cm)   | 13    |
| Packed unit / Weight      | 0.03 | Bare Unit / Weight (kg)   | 0.028 |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|

📅 18/06/2026 23:52:41 – For information purposes only, can be subject to change