

DEBAKEY FORCEPS 2.5 MM 15 CM
ATRAUMATIC MORS

Ref. IP27025



DEBAKEY FORCEPS 2.5 MM 15 CM ATRAUMATIC MORS

Ref. [IP27025](#)



Allows stitches or dissections to be performed

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104715226	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.6
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	17.7
Packed unit / Weight	0.035	Bare Unit / Weight (kg)	0.035

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---