

## DOYEN NEEDLE HOLDER 14 CM

Ref. IP30414



**DOYEN NEEDLE HOLDER 14 CM****Ref. IP30414****Features**

| CE Marking        | CE             | Medical device | Class I  |
|-------------------|----------------|----------------|----------|
| Bar code          | 08961104719040 | HS Code        | 90189084 |
| Warranty (months) | 24             | Country        | Pakistan |

**DIMENSIONS AND WEIGHT**

|                           |       |                           |       |
|---------------------------|-------|---------------------------|-------|
| Packed unit / Width (cm)  | nc    | Packed unit / Height (cm) | nc    |
| Packed unit / Length (cm) | nc    | Bare Unit / Length (cm)   | 14    |
| Packed unit / Weight      | 0.036 | Bare Unit / Weight (kg)   | 0.036 |

**SPECIFIC STORAGE CONDITIONS**

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

**OTHER FEATURES**

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|