

DOYEN NEEDLE HOLDER 15 CM

Ref. IP30415



DOYEN NEEDLE HOLDER 15 CM

Ref. IP30415



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104719057	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15
Packed unit / Weight	0.039	Bare Unit / Weight (kg)	0.039

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---