

OLSEN-HÉGAR NEEDLE HOLDER 14 CM

Ref. IP32414



OLSEN-HÉGAR NEEDLE HOLDER 14 CM

Ref. IP32414



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104719071	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	14
Packed unit / Weight	0.034	Bare Unit / Weight (kg)	0.034

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---