

CRAFOORD FORCEPS 24 CM

Ref. IP34600



Features

| | | | |
|------------------------------------|----------------|-------------------------|----------|
| CE Marking | CE | Medical device | Class Ir |
| Bar code | 08961104716742 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 8 | Bare Unit / Width (cm) | 8.6 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 0.4 |
| Packed unit / Length (cm) | 28 | Bare Unit / Length (cm) | 28 |
| Packed unit / Weight | 0.076 | Bare Unit / Weight (kg) | 0.07 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

📅 19/06/2026 02:18:42 – For information purposes only, can be subject to change