

WHITE TONSIL FORCEPS 23 CM

Ref. IP44223



Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104720466	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	23.2
Packed unit / Weight	0.066	Bare Unit / Weight (kg)	0.066

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---