

MEMBRANE PUNCH 19 CM

Ref. IP44401



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104721043	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	5	Bare Unit / Width (cm)	5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.1
Packed unit / Length (cm)	19	Bare Unit / Length (cm)	19
Packed unit / Weight	0.015	Bare Unit / Weight (kg)	0.015

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---