

## HUNTER SPLINT FORCEPS STRAIGHT 10.5 CM

Ref. IP45510



# HUNTER SPLINT FORCEPS STRAIGHT 10.5 CM

Ref. IP45510



## Features

|                   |                |                |          |
|-------------------|----------------|----------------|----------|
| CE Marking        | NONE           | Medical device | No       |
| Bar code          | 08961104714342 | HS Code        | 90189084 |
| Warranty (months) | 24             | Country        | Pakistan |

### DIMENSIONS AND WEIGHT

|                           |       |                           |       |
|---------------------------|-------|---------------------------|-------|
| Packed unit / Width (cm)  | nc    | Packed unit / Height (cm) | nc    |
| Packed unit / Length (cm) | nc    | Packed unit / Weight      | 0.012 |
| Bare Unit / Weight (kg)   | 0.012 |                           |       |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES

|                        |   |  |  |
|------------------------|---|--|--|
| Packed unit / Quantity | 1 |  |  |
|------------------------|---|--|--|