

BALLENGER SPONGE FORCEPS STRAIGHT 18 CM (STANDARD EYE)

Ref. IP45618



BALLENGER SPONGE FORCEPS STRAIGHT 18 CM (STANDARD EYE)

Ref. IP45618



Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104720114	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7.2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	17.7
Packed unit / Weight	0.046	Bare Unit / Weight (kg)	0.046
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		