

GUYON FORCEPS 23 CM

Ref. IP46823



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104732025	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Bare Unit / Width (cm)	7.6
Packed unit / Height (cm)	3	Bare Unit / Height (cm)	3
Packed unit / Length (cm)	22	Bare Unit / Length (cm)	22.7
Packed unit / Weight	0.074	Bare Unit / Weight (kg)	0.074
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

18/06/2026 21:53:39 – For information purposes only, can be subject to change