

DUJARDIN EYE BALL FIXATION FORCEPS, 11.5CM

Ref. IP52100



Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104714687	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	20
Packed unit / Weight	0.013	Bare Unit / Weight (kg)	0.013
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 19/06/2026 02:51:38 – For information purposes only, can be subject to change