

## GUNNAR DRAWING FORCEPS 18.5 CM

Ref. IP54818



CE

DISPOSITIF MEDICAL  
CLASSE  
**Ir**  
MEDICAL DEVICE

GARANTIE  
**2**  
2 YEARS  
GUARANTEE

STAINLESS  
**INOX**  
STEEL

19/06/2026 01:20:19 – For information purposes only, can be subject to change

# GUNNAR DRAWING FORCEPS 18.5 CM

Ref. IP54818



## Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104720732	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.8
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	18.4
Packed unit / Weight	0.098	Bare Unit / Weight (kg)	0.098
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 19/06/2026 01:20:19 – For information purposes only, can be subject to change