

## BAR UMBILICAL CLAMP CUTTING FORCEPS, SEMI-SHARPED

Ref. IP56401



# BAR UMBILICAL CLAMP CUTTING FORCEPS, SEMI-SHARPED

Ref. IP56401



## Features

|                             |                |                         |          |
|-----------------------------|----------------|-------------------------|----------|
| CE Marking                  | CE             | Medical device          | Class Ir |
| Bar code                    | 08961104737969 | HS Code                 | 90189084 |
| Warranty (months)           | 24             | Country                 | Pakistan |
| DIMENSIONS AND WEIGHT       |                |                         |          |
| Packed unit / Width (cm)    | nc             | Bare Unit / Width (cm)  | 3.6      |
| Packed unit / Height (cm)   | nc             | Bare Unit / Height (cm) | 2.5      |
| Packed unit / Length (cm)   | nc             | Bare Unit / Length (cm) | 19.2     |
| Packed unit / Weight        | 0.112          | Bare Unit / Weight (kg) | 0.112    |
| SPECIFIC STORAGE CONDITIONS |                |                         |          |
| Minimum temperature         | NA             | Max temperature         | NA       |
| Minimum humidity            | NA             | Maximum humidity        | NA       |
| OTHER FEATURES              |                |                         |          |
| Packed unit / Quantity      | 1              |                         |          |