

RAZOR, HANDLE STAINLESS STEEL

Ref. IR04100



RAZOR, HANDLE STAINLESS STEEL

Ref. IR04100



Features

CE Marking	NONE	Medical device	No
Bar code	08961104715738	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15.6
Packed unit / Weight	0.053	Bare Unit / Weight (kg)	0.053
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		