

STAINLESS STEEL RULER 30 CM

Ref. IR10030



STAINLESS STEEL RULER 30 CM

Ref. IR10030



Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104715769	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	3	Bare Unit / Width (cm)	3.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.1
Packed unit / Length (cm)	33	Bare Unit / Length (cm)	33.2
Packed unit / Weight	0.017	Bare Unit / Weight (kg)	0.017

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---