

HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°1/2

Ref. IS01401

Exemple de bougies Hégar double, inox.
15 tailles disponibles, nous consulter.



HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°1/2

Ref. IS01401



Description

Hégar Double Ended Dilator Stainless Steel

15 Available sizes, contact us

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104736832	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	4.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19
Packed unit / Weight	0.007	Bare Unit / Weight (kg)	0.007

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---