

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°23

Ref. IS01623



Exemple de bougies Hégar simple, inox.
30 tailles disponibles, nous consulter



HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°23

Ref. IS01623



Description

Hégar Single Ended Dilator Stainless Steel br>
30 Sizes available, contact us

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104736726	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19
Packed unit / Weight	0.125	Bare Unit / Weight (kg)	0.125

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---