

HÉGAR SINGLE ENDED DILATOR S. S. 19 CM SET OF 15 N°4 TO 18



Ref. IS01631

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104736818	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	10.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6.2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	21.5
Packed unit / Weight	0.787	Bare Unit / Weight (kg)	0.787
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		