

**METAL FEMALE CATHETER 15 CM X 4 MM**Ref. IS03400**Features**

|                             |                |                         |          |
|-----------------------------|----------------|-------------------------|----------|
| CE Marking                  | CE             | Medical device          | Class I  |
| Bar code                    | 08961104732469 | HS Code                 | 90189084 |
| Warranty (months)           | 24             | Country                 | Pakistan |
| DIMENSIONS AND WEIGHT       |                |                         |          |
| Packed unit / Width (cm)    | nc             | Bare Unit / Width (cm)  | 0.5      |
| Packed unit / Height (cm)   | nc             | Bare Unit / Height (cm) | 0.5      |
| Packed unit / Length (cm)   | nc             | Bare Unit / Length (cm) | 16       |
| Packed unit / Weight        | 0.006          | Bare Unit / Weight (kg) | 0.006    |
| SPECIFIC STORAGE CONDITIONS |                |                         |          |
| Minimum temperature         | NA             | Max temperature         | NA       |
| Minimum humidity            | NA             | Maximum humidity        | NA       |
| OTHER FEATURES              |                |                         |          |
| Packed unit / Quantity      | 1              |                         |          |