

**PINARD OBSTETRICAL STETHOSCOPE,
ALUMINIUM**

Ref. IS07601



19/06/2026 08:20:02 – For information purposes only, can be subject to change

PINARD OBSTETRICAL STETHOSCOPE, ALUMINIUM

Ref. IS07601



Stethoscope foetal Pinard used to listen to foetal heart sounds

Description

Small conical instrument, of which the top allows to the practitioner to place the ear, while the flat and circular basis is placed on the abdomen of the pregnant woman.
Designed for frequent and disinfection with hospital-grade products.

- Single chestpiece
- Lightweight
- Unbreakable aluminum
- Easy to handle
- Easy to clean

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104715752	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	13	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	7	Bare Unit / Height (cm)	7
Packed unit / Length (cm)	23	Bare Unit / Length (cm)	14
Packed unit / Weight	0.034	Bare Unit / Weight (kg)	0.034

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 19/06/2026 08:20:02 – For information purposes only, can be subject to change