

PINARD STETHOSCOPE WOOD

Ref. IS07602



CE



19/06/2026 07:54:27 – For information purposes only, can be subject to change

PINARD STETHOSCOPE WOOD

Ref. IS07602



Description

Packaging:

Individual high-density polyethylene bag

Outer packaging:

Cardboard box (box of 10)

Storage and Handling:

The Pinard stethoscope (wooden, 15 cm) must be handled with care and stored in a clean, dry place, away from direct sunlight, excessive humidity, and extreme temperatures. To preserve its quality, shape, and the smooth finish of the wood, must be stored separately and not placed under heavy objects or stacked with other instruments, as this may cause cracks, warping, or surface damage.

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104715745	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15.1
Packed unit / Weight	0.054	Bare Unit / Weight (kg)	0.052

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 19/06/2026 07:54:27 – For information purposes only, can be subject to change