

## EAR SYRINGE 110 ML

Ref. IS17110



### Features

|                   |                |                |          |
|-------------------|----------------|----------------|----------|
| CE Marking        | CE             | Medical device | Class I  |
| Bar code          | 08961104729155 | HS Code        | 90189084 |
| Warranty (months) | 24             | Country        | Pakistan |

#### DIMENSIONS AND WEIGHT

|                           |      |                         |       |
|---------------------------|------|-------------------------|-------|
| Packed unit / Width (cm)  | nc   | Bare Unit / Width (cm)  | 6     |
| Packed unit / Height (cm) | nc   | Bare Unit / Height (cm) | 5.5   |
| Packed unit / Length (cm) | nc   | Bare Unit / Length (cm) | 23    |
| Packed unit / Weight      | 0.28 | Bare Unit / Weight (kg) | 0.279 |

#### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

#### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|