

TAP WITH "T" HANDLE DIAM 3.5 MM S/S

Ref. IT00001



Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104733237	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

19/06/2026 07:07:51 – For information purposes only, can be subject to change