

LOUTE LOOP TWISTER

Ref. IT10100



LOUTE LOOP TWISTER

Ref. IT10100



Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104732773	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Bare Unit / Width (cm)	7.8
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	4.3
Packed unit / Length (cm)	26	Bare Unit / Length (cm)	26
Packed unit / Weight	0.378	Bare Unit / Weight (kg)	0.378
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		