

OCHSNER TROCART 14 FRENCH GAUGE (4.66 MM)

Ref. IT11614



Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104725256	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	7	Packed unit / Height (cm)	1
Packed unit / Length (cm)	21	Packed unit / Weight	0.028
Bare Unit / Weight (kg)	0.028		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--