

STAINLESS STEEL BED PAN SLIPPER TYPE 42X29 CM WITHOUT LID

Ref. AI14832



📅 19/06/2026 08:22:38 – For information purposes only, can be subject to change

STAINLESS STEEL BED PAN SLIPPER TYPE 42X29 CM WITHOUT LID

Ref. AI14832



Description

- Maximum weight supported: 90 kg

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104739673	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	30	Bare Unit / Width (cm)	29
Packed unit / Height (cm)	17	Packed unit / Length (cm)	56
Bare Unit / Length (cm)	42	Packed unit / Weight	1.08
Bare Unit / Weight (kg)	1.481		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---