

LID FOR BEDPAN (STAINLESS STEEL)

Ref. AI14833



19/06/2026 04:24:01 – For information purposes only, can be subject to change

LID FOR BEDPAN (STAINLESS STEEL)

Ref. AI14833



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104739666	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	31	Packed unit / Height (cm)	2
Packed unit / Length (cm)	31	Packed unit / Weight	0.416
Bare Unit / Weight (kg)	0.416		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		