

S/S TRAY WITH COVER 280 X 180 X 30 MM

Ref. AI30214



S/S TRAY WITH COVER 280 X 180 X 30 MM

Ref. AI30214



Ideal for storing small surgical instruments

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104739871	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	21	Bare Unit / Width (cm)	21.5
Packed unit / Height (cm)	5	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	30	Bare Unit / Length (cm)	30
Packed unit / Weight	0.95	Bare Unit / Weight (kg)	0.95

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---