

STAINLESS STEEL INSTRUMENT TRAY 46.5 X 31 X 5 CM

Ref. AI46315



Ideal for storing small surgical instruments

STAINLESS STEEL INSTRUMENT TRAY 46.5 X 31 X 5 CM

Ref. AI46315



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104739857	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	33	Packed unit / Height (cm)	6
Packed unit / Length (cm)	48	Packed unit / Weight	1

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---