

STAINLESS STEEL INSTRUMENT TRAY 52 X 34 X 5 CM

Ref. AI46317



18/06/2026 22:37:53 – For information purposes only, can be subject to change

STAINLESS STEEL INSTRUMENT TRAY 52 X 34 X 5 CM

Ref. AI46317



Ideal for storing small surgical instruments

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104749030	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	5	Bare Unit / Width (cm)	34
Packed unit / Height (cm)	33	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	47	Bare Unit / Length (cm)	52
Packed unit / Weight	1		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---