

LEXER GOUGE 25CM 4MM

Ref. IG01131



Features

CE Marking	CE	Medical device	Class Ir
Bar code	08961104733343	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	2.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	25
Packed unit / Weight	0.109	Bare Unit / Weight (kg)	0.109
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		