

# **MINI-HOHMANN BONE LEVER 17 CM X 6 MM**

**Ref. IL11818**



## **Features**

|                   |                |                |          |
|-------------------|----------------|----------------|----------|
| CE Marking        | CE             | Medical device | Class I  |
| Bar code          | 08961104735231 | HS Code        | 90189084 |
| Warranty (months) | 24             | Country        | Pakistan |

### **DIMENSIONS AND WEIGHT**

|                           |       |                         |       |
|---------------------------|-------|-------------------------|-------|
| Packed unit / Width (cm)  | 1     | Bare Unit / Width (cm)  | 1.4   |
| Packed unit / Height (cm) | nc    | Bare Unit / Height (cm) | 0.3   |
| Packed unit / Length (cm) | 17    | Bare Unit / Length (cm) | 17.2  |
| Packed unit / Weight      | 0.025 | Bare Unit / Weight (kg) | 0.025 |

### **SPECIFIC STORAGE CONDITIONS**

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### **OTHER FEATURES**

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|