

**KOCHER FORCEPS STRAIGHT 1X2 TH 18  
CM ( COLLIN LOCK)**

Ref. IP02418



CE



19/06/2026 07:16:26 – For information purposes only, can be subject to change

# KOCHER FORCEPS STRAIGHT 1X2 TH 18 CM ( COLLIN LOCK)



Ref. IP02418

## Features

CE Marking	CE	Medical device	Class I
Bar code	08961104716612	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.037
Bare Unit / Weight (kg)	0.037		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

19/06/2026 07:16:26 – For information purposes only, can be subject to change