

# DEBAKEY FORCEPS 1,5 MM X 12 CM ATRAUMATIC JAWS



Ref. IP27112

**Allows stitches or dissections to be performed**

## Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations.

Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016).

- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" leaflet

## Features

| CE Marking        | CE             | Medical device | Class I  |
|-------------------|----------------|----------------|----------|
| Bar code          | 08961104715165 | HS Code        | 90189084 |
| Warranty (months) | 24             | Country        | Pakistan |

### DIMENSIONS AND WEIGHT

|                           |       |                           |       |
|---------------------------|-------|---------------------------|-------|
| Packed unit / Width (cm)  | nc    | Packed unit / Height (cm) | nc    |
| Packed unit / Length (cm) | nc    | Packed unit / Weight      | 0.029 |
| Bare Unit / Weight (kg)   | 0.029 |                           |       |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES

|                        |   |  |  |
|------------------------|---|--|--|
| Packed unit / Quantity | 1 |  |  |
|------------------------|---|--|--|