



CARDIOMATE PC ECG SPENGLER

Ref. 260039



CE
0123



19/06/2026 09:05:46 – For information purposes only, can be subject to change

CARDIOMATE PC ECG SPENGLER

Ref. 260039



A SPENGLER PC-ECG with XXL performance in the palm of your hand

Description

Digital ECG, the Cardiomate PC allows the real-time display of a 12-lead ECG on the screen of a computer and its printing on an A4 printer.

Compatible with a Windows system, the Cardiomate PC connects via a simple USB cable. With the same performance as the Cardiomate range and with its compact and lightweight case, the Cardiomate PC is a state-of-the-art measuring instrument in the office or on the go.

The software allows the visualization, recording, analysis and storage of ECGs in a clear and intuitive way.

VERY HIGH ACCURACY OF

MEASUREMENTS- Visualization of the 12 leads in real time on PC

- Sampling speed 32,000Hz

- 24-bit converter

- Automatic interpretation and calculation of intervals

- Pacemaker detection

EASY TO USE

- Connection to the PC by USB

- Self-powered- Simple

- and intuitive software

- Software compatible with other ECG Cardiomate

COMPACT AND LEGER

- Dimensions: 68x126x24mm

- Weight: 104g

CARDIOMATE PC ECG SPENGLER

Ref. 260039



Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446042059	HS Code	9018110000
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	12.6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	6.7
Packed unit / Weight	nc	Bare Unit / Weight (kg)	0.104

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---