


PULSE STETHOSCOPE / RED

Ref. 507521



CE



📅 19/06/2026 07:39:37 – For information purposes only, can be subject to change

PULSE STETHOSCOPE / RED

Ref. 507521



Description

For daily blood pressure monitoring. Powder-coated or chromium-plated (black model) zinc alloy single head chestpiece. Epoxy resin, high resolution acoustic membrane. Screw-on flexible tips to match the PVC tubing. Angled binaura with built-in spring.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009274	HS Code	90181990
Warranty (months)	36		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.22
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Azure	Packed unit / Quantity	1

📅 19/06/2026 07:39:37 – For information purposes only, can be subject to change