

## MAXI +3 WHEEL STAND ADULT CUFF (M)

Ref. 522705



CE  
0123

DISPOSITIF MEDICAL  
CLASSE  
**I**m  
MEDICAL DEVICE

GARANTIE  
2 ANS  
**2**  
2 YEARS  
GUARANTEE

📅 19/06/2026 03:18:24 – For information purposes only, can be subject to change

## MAXI +3 WHEEL STAND ADULT CUFF (M)

Ref. 522705



**SPENGLER technology in a blood pressure monitor with a dial with a diameter of 147mm.**

### Description

#### HIGH MEASUREMENT

ACCURACY• Double tubing for linear pressure rise

#### COMFORT• Ambidextrous•

Remotely readable measurement•

2.5-meter spiral tubing• Large

volume ergonomic pear for quick inflation•

Storage basket behind the dial• Comes

with cuff wheeled base

• Comes with armband size M, nylon

#### OPTION

• 3 sizes of cuffs (S, M, L) available

# MAXI +3 WHEEL STAND ADULT CUFF (M)

Ref. 522705



## Features

|                                    |               |                           |                               |
|------------------------------------|---------------|---------------------------|-------------------------------|
| CE Marking                         | CE0123        | Medical device            | Class Im (measuring function) |
| Bar code                           | 3700446028152 | HS Code                   | 90189010                      |
| Warranty (months)                  | 24            | Country                   | China                         |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |                               |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc                            |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 4.16                          |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |                               |
| Minimum temperature                | NA            | Max temperature           | NA                            |
| Minimum humidity                   | NA            | Maximum humidity          | NA                            |
| <b>OTHER FEATURES</b>              |               |                           |                               |
| Packed unit / Quantity             | 1             |                           |                               |

19/06/2026 03:18:24 – For information purposes only, can be subject to change