

MAXI +3 WALL MOUNTED WITH SET OF 3 CUFFS (S, M, L)

Ref. 522707



CE
0123



MAXI +3 WALL MOUNTED WITH SET OF 3 CUFFS (S, M, L)



Ref. 522707

SPENGLER technology in a blood pressure monitor with a dial with a diameter of 147mm.

Description

HIGH MEASUREMENT

ACCURACY• Double tubing for linear pressure rise

COMFORT• Ambidextrous•

Remotely readable measurement•

2.5-meter spiral tubing• Large volume ergonomic pear for quick inflation•

Storage basket behind the dial•

Adjustable wall mount

ARMBAND• Supplied with kit of 3 cuffs (S, M, L), nylon

Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3700446031152	HS Code	90189010
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.154

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---